Colour Reference:

BLACK

P 689 C

P 287 U

PLEASE NOTE:

***This document to be
exported for proof and
printing as pairs because
the actual page size is that
of two single pages
i.e. 150 x 105mm***

Colour Reference:

BLACK

P 689 C

P 287 U

INSTRUCTIONS FOR USE

1

Phase-pure, porous silicated calcium phosphate granules for use as a synthetic bone graft substitute.

Sterile - Single use - Store in a cool, dry place.

Description:

- 0.8 wt% Silicated calcium phosphate
- Granule size 1 – 2 mm
- Total porosity 80% ± 2.5%
- Aqueous carrier

Indications:

- Bone graft substitutes are intended to be used in place of corticocancellous, or cancellous allograft or autograft bone.
- The mechanical environment for such uses experience either low load requirements or compression.
- Typical surgical applications for bone graft substitutes are:-
 - Small void filling, e.g. after removal of a small bone tumour or following bone fracture reduction or in osteotomies.
 - Spinal fusion, where a cage or screw fixation device is used to relieve the graft site from physiological loads.
- It is not intended to be used in place of cortical strut allograft bone where high tensile, torsion and/or bending strength are required. The products are used by surgeons in place of allograft bone (bone from humans stored in bone banks).

Directions for Use:

- Use in place of/or added to morsellised autograft bone.
- Use in place of/or added to morsellised allograft bone.
- Actifuse ABX is designed to be used alone. However, it may be mixed with sterile saline, autologous blood or bone marrow aspirate but this may have an effect on handling.
- Actifuse ABX implanted alone or when mixed is used in the same manner as morsellised allograft or autograft.

Risks and warnings:

- The risks with Actifuse ABX (e.g. infection, graft loosening, non-union, delayed union, repeat surgery) are the same as those for allograft bone. However, with Actifuse ABX there is no risk of disease transmission.

Colour Reference:

BLACK

P 689 C

P 287 U

Colour Reference:

BLACK

P 689 C

P 287 U

- Actifuse ABX should not be used where it could be subject to tension, torsion, compression, shear or bending. A conventional implant (e.g., screw, rod) can protect the graft from such loading actions.
- Actifuse ABX should not be used in volumetrically unconstrained sites (so the graft material cannot move or escape).
- Do not overfill or attempt to pressurize the bony defect site, as this may lead to extrusion of the product beyond the site of its intended application and damage to the surrounding tissues, or may lead to fat embolization or embolization of the product into the blood stream.
- Additional mixing containers and stirrers must be sterilized before use.
- Actifuse should be used with caution in pediatric cases harboring large Juvenile Bone Cysts. Isolated cases of transient postoperative fever and inflammatory reaction, in the absence of infection have been reported from clinical experience with extremely large cysts. This side effect had no negative impact on the therapeutic outcome.
- The safety and effectiveness of the Actifuse Bone Graft Substitutes has not been studied in children under 8 and thus is not known.

Sterility:

- Sterilised by gamma irradiation.
- Do not resterilise.
- If the packaging or seal are damaged do not use. Discard.
- Discard unused product. Re-use of single-use devices creates a potential risk of patient or user infections.

Contraindications:

- Impaction grafting for failed total hip or knee arthroplasty
- Direct loading of graft material, i.e. in the absence of conventional implants, such as screws and rods
- Infection
- Inability to cover or deliberate non-coverage of graft site using soft tissue
- Avascular or compromised vascular network sites
- Patient metabolism may compromise bony regeneration
- Avoid use in patients where in the surgeon's opinion patient lifestyle, compliance and/or physical attributes would compromise clinical outcome
- Medication which could slow bone healing

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Baxter Healthcare Corporation
Deerfield, IL 60015 USA

Baxter and ACTIFUSE are registered trademarks of Baxter International Inc., or its subsidiaries.

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Revision number 0723264, 2019-10-01

Baxter

Release Status: Issued and Effective

Colour Reference:

BLACK

P 689 C

P 287 U

Colour Reference:

BLACK

P 689 C

P 287 U

Store between 5°C-32°C
 Zwischen 5 °C und 32 °C lagern.
 Φυλάσσεται σε θερμοκρασία 5°C-32°C
 Almacenar entre 5°C y 32°C
 Oppbevares mellom 5 °C og 32 °C
 Armazenar entre 5°C-32°C
 Säilytettävä 5–32 °C:n lämpötilassa
 Conservare tra 5 °C e 32 °C
 A conserver entre 5 °C et 32 °C
 Opbevares mellem 5°C-32°C
 Förvara vid en temperatur på 5°C-32°C
 Bewaren tussen 5°C en 32°C
 32°C ו 5°C אחסן בטמפרטורה שבין 5°C-32°C
 Przechowywać w temperaturze 5°C-32°C
 Skladujte při teplotě 5 °C–32 °C
 Да се съхранява при температура 5°C-32°C
 Tárolás 5 °C–32 °C hőmérsékleten
 Hoida temperatuuril 5°C-32°C
 Laikyti 5–32 °C temp.
 A se depozita la temperaturi situate în intervalul 5°C-32°C
 Skladujte pri teplote 5 °C-32 °C

Not made with rubber natural latex
 Nicht aus Naturlatex hergestellt
 Κατασκευασμένο χωρίς φυσικό ελαστικό κόμμι (λάτεξ)
 No contiene látex de caucho natural
 Ikke produsert av naturgummilateks
 Não contém látex de borracha natural
 Ei sisällä luonnonkumilateksia
 Non contiene lattice di gomma naturale
 Ne contient pas de latex naturel
 Ikke fremstillet med latex naturgummi
 Innehåller inte naturgummi (latex)
 Bij de productie is geen natuurlijk rubber (latex) gebruikt
 אינו עשוי מגומי לטקס טבעי
 Nie zawiera kauczuku naturalnego (lateksu)
 Neobsahuje přírodní kaučukový latex
 Не е изработено от латекс от естествен каучук
 Nem természetes latex gumbóli készült
 Ei sisalda looduslikku lateksit
 Pagaminta be natūralios latekso gumos
 Nu este fabricat cu latex din cauciuc natural
 Neobsahuje prírodný kaučukový latex.

Colour Reference:

BLACK

P 689 C

P 287 U

Release Status:Issued and Effective

TcU ELECTRONIC SIGNATURE REPORT

REVISION INFORMATION				
Item ID: EL880723264		Revision ID: 1		
Item Name: ACTIFUSE ABX ROW IFU		Release Date: 17-Dec-2019		
Description: ACTIFUSE ABX ROW IFU				
CHANGE INFORMATION				
CN/CR Number (if applicable): N/A				
Description of Change (This field will be blank if required data is not available): New artwork for ABX IFU				
Reason for Change (This field will be blank if required data is not available): CC 1570887				
APPROVALS & SIGNATURES for Document Release				
Name	Role	Workflow Step	Date of Signature	Decision Taken
Gorasia, Priya Dinesh	SME	Initiate Review	10-Dec-2019	Approved
Wilkins, Rebecca	SME	Document Review - SME & Quality	10-Dec-2019	Approved
Mason, Paul John	Quality	Document Review - SME & Quality	16-Dec-2019	Approved
Kumar, Jaishree	Change Specialist 3	Release Document(s)	17-Dec-2019	Approved
Kumar, Jaishree	Change Specialist 3	Set Effectivity	17-Dec-2019	Approved

Release Status: Issued and Effective

Baxter Confidential - Internal Use Only