BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: BE-30-03-130	Date: 06FEB19	Proofread No.: P03
Designer: Delcourt Perrine	Page: 1 of 108	
Colour Reference:	BLACK	

HEMOPATCH

Sealing Hemostat

Resorbable Collagen-Based
Sealing Hemostat

Baxter

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: BE-30-03-130

Date: 06FEB19

Proofread No.: P03

Designer: Delcourt Perrine

Page: 2 of 108

Colour Reference:

BLACK

Contents

01	eng	HEMOPATCH SEALING HEMOSTAT
05	bul	ЗАПЕЧАТВАЩ ХЕМОСТАТИК HEMOPATCH
09	cze	TĚSNÍCÍ HEMOSTAT HEMOPATCH
13	dan	HEMOPATCH FORSEGLENDE HÆMOSTATIKUM
17	dut	HEMOPATCH SEALING HAEMOSTATICUM
21	est	HEMOPATCH HERMETISEERIV HEMOSTAATIKUM
25	fin	HEMOPATCH TIIVISTÄVÄ HEMOSTAATTI
29	fre	PATCH HÉMOSTATIQUE D'ÉTANCHÉITÉ HEMOPATCH
33	ger	HEMOPATCH HÄMOSTATIKUM MIT VERSIEGELNDEN EIGENSCHAFTEN
37	gre	ΑΙΜΟΣΤΑΤΙΚΟ ΣΤΕΓΑΝΟΠΟΙΗΣΗΣ HEMOPATCH
41	hrv	HEMOPATCH HEMOSTATIČKO LJEPILO
45	hun	HEMOPATCH LEZÁRÓ HAEMOSTATICUM
49	ice	HEMOPATCH BLÆÐINGAHEMJANDI SÁRAPÚÐI
53	ita	HEMOPATCH SIGILLANTE EMOSTATICO
57	lav	HEMOPATCH HERMETIZĒJOŠS HEMOSTATISKS LĪDZEKLIS
61	lit	„HEMOPATCH“ HERMETINĖ HEMOSTAZINĖ PLĖVELĖ
65	nor	HEMOPATCH FORSEGLENDE HEMOSTATIKUM
69	pol	HEMOPATCH HEMOSTATYK USZCZELNIĄCY
73	por	HEMOPATCH SELANTE HEMOSTÁTICO
77	rum	HEMOPATCH PLASTURE HEMOSTATIC
81	rus	ГЕРМЕТИЗИРУЮЩЕЕ ГЕМОСТАТИЧЕСКОЕ ИЗДЕЛИЕ HEMOPATCH
85	slo	HEMOSTATIKUM NA ZASTAVENIE KRVÁCANIA HEMOPATCH
89	slv	LEPILNI HEMOSTATIK HEMOPATCH
93	spa	HEMOPATCH SELLADO HEMOSTÁTICO
97	swe	HEMOPATCH FÖRSEGLANDE HEMOSTAS
101	tur	HEMOPATCH SIZDIRMAZ HEMOSTAT

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: BE-30-03-130	Date: 06FEB19	Proofread No.: P03
Designer: Delcourt Perrine	Page: 3 of 108	
Colour Reference:	BLACK	

HEMOPATCH SEALING HEMOSTAT
**RESORBABLE COLLAGEN-BASED
 SEALING HEMOSTAT**

DESCRIPTION

HEMOPATCH Sealing Hemostat ("HEMOPATCH") consists of a soft, thin, pliable, flexible pad of collagen derived from bovine dermis, coated with NHS-PEG (pentaerythritol polyethylene glycol ether tetra-succinimidyl glutarate).

Due to its flexible structure, the application of HEMOPATCH to the site where hemostasis / sealing is desired is easily controlled. For differentiation, the non-coated side is marked with blue squares using a biocompatible colorant.

HEMOPATCH is supplied in the following 3 sizes:

- HEMOPATCH 27 x 27 mm
- HEMOPATCH 45 x 45 mm
- HEMOPATCH 45 x 90 mm

INTENDED USE

HEMOPATCH is an absorbable collagen pad intended for sealing and hemostasis.

INDICATION

HEMOPATCH is indicated as a hemostatic device and surgical sealant for procedures in which control of bleeding or leakage of other body fluids or air by conventional surgical techniques is either ineffective or impractical. HEMOPATCH may be used to close dural defects following traumatic injury, excision, retraction or shrinkage of the dura mater.

MECHANISM OF ACTION

The performance of HEMOPATCH involves the interaction of two components that are able to create rapid and lasting hemostasis by sealing off the bleeding surface.

The white, tissue-facing side of HEMOPATCH is covered with a thin layer of NHS-PEG. The NHS-PEG, when in contact with blood or other body fluids forms a hydrogel which enhances its adhering properties and seals the tissue surface.

In contact with blood, collagen induces aggregation of platelets. Platelets deposit in large numbers on the collagen structure, degranulate, and release coagulation factors that, together with plasma factors, enable the formation of fibrin. The structure of HEMOPATCH provides a three dimensional matrix that readily absorbs fluids, and provides additional mechanical strengthening of the clot.

Preclinical studies demonstrated resorption and replacement with host tissue in 6 – 8 weeks with little tissue reaction. The rate of resorption and replacement with host tissue can be longer as it is dependent on several factors, including the amount of product left in situ and the site of use. When used for dura replacement, the collagen pad of HEMOPATCH serves as a scaffold for tissue ingrowth and will be replaced over time (in a preclinical study this time occasionally extended past 16 weeks) by newly formed tissue.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: BE-30-03-130

Date: 06FEB19

Proofread No.: P03

Designer: Delcourt Perrine

Page: 4 of 108

Colour Reference:

BLACK

02 eng

CONTRAINDICATIONS

Do not compress HEMOPATCH into blood vessels or use intravascularly.

The device must not be used in patients with known hypersensitivity to bovine proteins or brilliant blue (FD&C Blue No. 1 (Blue 1)).

WARNINGS

HEMOPATCH is not intended to be used in pulsatile, severe bleedings.

The use of HEMOPATCH is not recommended in the presence of an active infection.

When used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, the brain and/or cranial nerves, care should be exercised to avoid overpacking (collagens may expand upon absorption of liquid), creating the potential for neural damage.

HEMOPATCH is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for hemostasis and sealing.

PRECAUTIONS

For single use only. Do not re-sterilize. Reuse of single-use devices creates a potential risk of infection.

Do not store without aluminum pouch.

Discard HEMOPATCH if the packaging is damaged.

Do not apply the surface marked with blue squares to the application site. Do not wet HEMOPATCH prior to application.

HEMOPATCH is self-adherent and is not intended to be sutured.

Do not apply on a dry tissue surface or lesion.

NHS-PEG only forms an adhering hydrogel when in contact with wound fluid such as blood or lymphatic. In the absence of such wound fluids, sodium bicarbonate solution¹ (concentration between 4.2% to 8.4%) can be used to moisten the tissue prior to application of HEMOPATCH.

Do not use HEMOPATCH on bone surfaces where adhesives, such as methylmethacrylate or other acrylic adhesives, will be required to attach a prosthetic device. Microfibrillar collagen has been reported to reduce the strength of methyl methacrylate adhesives used to attach prosthetic devices to bone surfaces.

The safety and performance of HEMOPATCH has not been established in children, pregnant and/or lactating women.

ADVERSE EVENTS

The following adverse events have been identified during the use of HEMOPATCH:

- Seroma formation
- Inflammatory reaction.

These reactions are reported voluntarily from a population of uncertain size and it is not always possible to reliably estimate their frequency or establish a causal relationship to device exposure.

Collagen-Based Agents Adverse Events:

Adverse events which may be possibly or probably related to the use of any collagens include:

- Seroma
- Hematoma
- Potentiation of infection
- Inflammatory reaction
- Foreign body reaction
- Adhesion formation
- Allergic reaction.

¹ the grade of sodium bicarbonate solution to be suitable for parenteral human use.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: BE-30-03-130

Date: 06FEB19

Proofread No.: P03

Designer: Delcourt Perrine

Page: 5 of 108

Colour Reference:

BLACK

eng 03

DRUG INTERACTIONS

Chlorine-releasing antiseptic agents as well as tannins and caustics should not be used in combination with collagen at the application site.

PREPARATION

HEMOPATCH comes ready to use in sterile packages and must be handled accordingly. Use only undamaged packages. Single use only. Do not re-sterilize.

Circulating Nurse:

- Open outer, non-sterile package (aluminum pouch) and,
- present the inner, sterile blister pack to the scrub nurse.

Scrub Nurse:

- Deliver the inner, sterile package into the sterile field, using aseptic technique. Once removed from the outer pouch HEMOPATCH can be used throughout surgery provided it remains dry and sterile.
- Open the inner, sterile package and present the device to the surgeon.
- HEMOPATCH should be maintained dry prior to application. Do not soak in solution prior to application (e.g. physiological saline solution, antiseptic agents).
- Dry gloves and instruments (e.g. forceps, scissors) should be used to handle HEMOPATCH.
- In situations in which HEMOPATCH has been removed for application from the inner pouch, but application is delayed, it can be used throughout surgery, provided it remains dry and sterile.

ADMINISTRATION

- **Do not apply intravascularly.**
- Select the appropriate size of the pad so that it overlaps the margins of the bleeding surface or tissue lesion by about 1 cm. The pad may be cut to the desired size and shape.
- Dry gloves and surgical instruments (forceps, scissors) should be used to handle, cut and apply HEMOPATCH.
- Multiple pads may be used for larger bleeding surfaces or tissue lesions. When using multiple pads, ensure a 1 cm overlap where the pads contact each other.

Maximum number of patches to be used

HEMOPATCH size	Adults	Children (≥10 Kg)
45 x 90 mm	7	1
45 x 45 mm	14	2
27 x 27 mm	28	4

Method of Application

- Apply the dry HEMOPATCH with the non-marked white surface contacting the bleeding or moist surface of the wound or tissue lesion. Improved adherence has been observed in preclinical studies when HEMOPATCH is in direct contact with wound fluid such as blood or lymphatic.
- Use a dry gauze or patty and hold in place with gentle, uniform pressure over the entire pad surface for 2 minutes.
- When applying HEMOPATCH, minimize contact with bloody or moist surgical instruments, gauzes or gloves due to the affinity of collagen to blood and body fluids.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: BE-30-03-130

Date: 06FEB19

Proofread No.: P03

Designer: Delcourt Perrine

Page: 6 of 108

Colour Reference:

BLACK

04 eng

- Gently remove the gauze or patty from the pad. Gentle irrigation may also help in removing the gauze or patty without dislodging HEMOPATCH from the surgical site.
- Leave HEMOPATCH in situ after hemostasis or sealing has been achieved.
- If hemostasis has not been satisfactorily achieved, HEMOPATCH can be removed, without tissue damage, up to 3 minutes from the time of the initial tissue contact. Use care when removing the pad. A new pad can then be reapplied according to the steps above.
- When using HEMOPATCH for sealing or for dura repair, removal of the pad, any time after application, can only be done at the surgeon's discretion.

SHELF LIFE AND STORAGE

Store in a dry place at 2° to 8°C.

You may remove the product from the refrigerator and store at room temperature (maximum 25°C) for up to six months within the indicated shelf life. Mark the date the product is taken out of the refrigerator (year/month/day) as shown on the symbol below that is printed on the shelf carton. Leave pouches in the shelf carton or mark the same date on the pouches as well.

YYYY/MM/DD

After storage at room temperature, HEMOPATCH should not be refrigerated again.

Do not re-sterilize.

Expiry date: See imprint on outer package.

Do not use after the expiry date printed on the label.

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Baxter and Hemopatch are registered trademarks of Baxter International Inc.