
# Peri-Strips Dry®

With Veritas® Collagen Matrix • Staple Line Reinforcement  
 Avec matrice de collagène Veritas® • Renforcement de ligne d'agrafes  
 Mit Veritas® Collagen-Matrix • Klammernaht-Verstärkung  
 Con matrice di collagene Veritas® • Rinforzo della linea di sutura  
 Con matriz de colágeno Veritas® • Refuerzo de la línea de grapado  
 Met Veritas® Collagen Matrix • Nietierlijnversterking  
 Med Veritas® Kollagen -matrix • Hæftesømsforstærkning  
 Med Veritas® Kollagenmatrix • Agraftlinjeförstärkning  
 Med Veritas® kollagen-matrise • stiftelinjeförsterkning  
 Veritas® Kollajen Matriksli • Zimba Hattı Kuvvetlendirici  
 με μήτρα κολλαγόνου Veritas® • Υλικό ενίσχυσης γραμμής συνδετήρων συρραπτικού  
 Veritas® Kollageenimatriisilla • Hakasrivin vahvike  
 Com Matriz de Colagénio Veritas® • Reforço da Linha de Agrafo  
 Cu matrice de colagen Veritas® • Intărirea liniei de sutură  
 S kolagenovou maticí Veritas® • Zpevnění linie svorek  
 Veritas® kollagénmátrixszal • Kapocssor-megerősítő  
 Z kolagenową macierzą Veritas® • Wzmocnienie linii szwu mechanicznego

<b>INSTRUCTIONS FOR USE</b> .....	<b>4</b>
<b>MODE D'EMPLOI</b> .....	<b>7</b>
<b>GEBRAUCHSANLEITUNG</b> .....	<b>10</b>
<b>ISTRUZIONI PER L'USO</b> .....	<b>14</b>
<b>INSTRUCCIONES DE USO</b> .....	<b>17</b>
<b>GEBRUIKSAANWIJZING</b> .....	<b>20</b>
<b>BRUGSVEJLEDNING</b> .....	<b>23</b>
<b>BRUKSANVISNING</b> .....	<b>26</b>
<b>BRUKSANVISNING</b> .....	<b>29</b>
<b>KULLANMA TALIMATI</b> .....	<b>32</b>
<b>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ</b> .....	<b>35</b>
<b>ΚΑΥΤΤÖÖHJEET</b> .....	<b>39</b>
<b>INSTRUÇÕES DE UTILIZAÇÃO</b> .....	<b>42</b>
<b>INSTRUCȚIUNI DE FOLOSIRE</b> .....	<b>45</b>
<b>NÁVOD K POUŽITÍ</b> .....	<b>48</b>
<b>HASZNÁLATI ÚTMUTATÓ</b> .....	<b>51</b>
<b>INSTRUKCJA UŻYCIA</b> .....	<b>54</b>

## Peri-Strips Dry® with Veritas® Collagen Matrix Components

## Peri-Strips Dry® with Veritas® Collagen Matrix Figures 1-7


FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4


FIGURE 5


FIGURE 6


FIGURE 7

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## SYMBOL DEFINITIONS:

-  Store at a controlled room temperature.
-  Keep away from heat. Do not use if heat indicator is red.
-  Do not reuse
-  Consult *Instructions for Use*
- STERILE** Sterile
- STERILE EO** Sterilized using ethylene oxide
-  Do not use if the product sterile barrier or its packaging is compromised.
- STERILE SH** Sterilized using steam heat
- NaOH** This product is treated with sodium hydroxide.
- BOVINE** This product is derived from USDA inspected cattle.
- MADE IN THE U.S.A.** Made in the USA
- Rx Only** CAUTION: Federal (USA) law restricts this device to sale by, or on the order of, a physician.
- REF** Catalog number
-  Use by date
- LOT** Lot number
-  Compatible stapler models
- REMOVE BEFORE USE** Remove before use
-  Manufacturer
- EC REP** Authorized Representative in the European Community
- PN** SSI part number
- TN** SSI tracking number
- MODEL** Model number
- ID** SSI internal code
-  Caution: Consult instruction for use for warning and precaution information.

## DESCRIPTION:

Peri-Strips Dry® with Veritas® Collagen Matrix (PSD-V) is prepared from dehydrated bovine pericardium procured from cattle under 30 months of age in the United States.

One (1) tube of PSD Gel (Gel) is provided for every two (2) pouches of PSD-V. The Gel is used to create a temporary bond between the PSD-V buttress and the surgical stapler jaws until the stapler is positioned and fired. Each PSD-V loading unit and each Gel tube is packaged sterile in a separate pouch.

PSD-V utilizes animal tissue; patient must be informed prior to any procedure.

## INDICATIONS FOR USE:

PSD-V is intended for use as a prosthesis for the surgical repair of soft tissue deficiencies using surgical staplers when staple line reinforcement is needed.

PSD-V can be used for reinforcement of staple lines during bariatric surgical procedures.

PSD-V can be used for reinforcement of staple lines during gastric, small bowel, colon, and colorectal procedures.

## CONTRAINDICATIONS:

The use of PSD-V is contraindicated in patients with known sensitivity to bovine material.

## ADVERSE REACTIONS:

As with any surgical procedure, adverse reactions are possible and include but are not limited to: infection, rejection, erosion, and allergic reaction.

## WARNINGS:

Do not re-sterilize. Resterilization may cause changes to the tissue and negatively impact functionality of the device.

Discard all open unused components.

Do not use product if there is damage to the pouch or seals.

Ensure the staple line is completely covered with the buttress, or inadequate coverage after firing may result.

PSD-V is not designed, sold, or intended for use except as indicated; doing so may result in surgical complications.

Synovis products differ; substitution of one product for another product may reduce product performance.

Do not use PSD-V if the product has been exposed (1) to solutions above room temperature, (2) to chemicals, antibiotics, or other substances other than specifically addressed in the instructions, (3) or if heat indicator has been activated.

#### CAUTIONS:

Use care when removing loading unit components from the stapler to prevent buttress dislodgement.

Do not get the buttress wet before applying Gel, or the buttress may not adhere to the stapler properly.

Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the buttress may not adhere to the stapler properly.

Final tissue compression, including PSD-V, must meet the range specified by the stapler manufacturer; this is especially important if staple firings are overlapped. PSD-V increases the total thickness of the area stapled by 0.4 mm - 1.2 mm (0.016" - 0.048").

Follow *Instructions for Use* supplied by the stapler manufacturer. Do not use PSD-V contrary to the stapler manufacturer's instructions.

The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.

#### INSTRUCTIONS FOR USE:

##### LOADING STAPLER FOR ALL STAPLER MODELS

*Note: Each model of PSD-V has been designed specifically for the stapler models indicated on the label; verify that the correct model of PSD-V has been selected.*

*Note: Loading technique for PSD-V varies. Follow the appropriate technique as indicated below.*

1. Open the outer PSD-V pouch.
2. Aseptically remove the inner pouch and place in the sterile field.

3. Inspect the pouch. Do not use if the pouch is damaged or if the seals are not intact.
  4. Open the inner PSD-V pouch and remove the PSD-V loading unit by using atraumatic techniques.
  5. Inspect the Gel pouch. Do not use if the pouch is damaged or if the seals are not intact.
  6. Open the Gel pouch and aseptically place Gel in the sterile field.
  7. Break off the metal tip of the Gel tube.
  8. Attach the plastic Gel Tip to the Gel tube by pushing it down over the top of the Gel tube.
  9. Apply a continuous thin Gel bead onto each buttress strip (See Figure 1).
 

**Caution:** Do not get the buttress wet before applying Gel or the buttress may not adhere to the stapler properly.
  10. Ensure that the jaws of the stapler are clean and dry before inserting the PSD-V loading unit.
 

**Caution:** Final tissue compression, including PSD-V, must meet the range specified by the stapler manufacturer. This is especially important if stapler firings are overlapped. PSD-V increases the total thickness of the area stapled by 0.4mm - 1.2mm (0.016" - 0.048").
  11. Identify the anvil (ANV) and cartridge (CART) sides of the PSD-V loading unit.
 

**Caution:** The cartridge and anvil loading units differ; substitution of one unit for another may interfere with closing and firing of the stapler.
- Note: For ENDO GIA™ UNIVERSAL stapler models skip to step 19.**
- ##### LOADING STAPLER FOR ALL STAPLER MODELS EXCLUDING ENDO GIA™ UNIVERSAL MODELS
12. Position the open stapler onto the stapler loading unit (See Figure 2).
  13. Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the buttress may not adhere to the stapler properly.
 

**Note:** The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.
  14. Close the stapler.
  15. Remove the plastic sheath leaving the foam spacer between the stapler jaws (See Figure 3).

16. Discard plastic sheath and remove any excess Gel from stapler.

**Note:** *PSD-V can be used immediately or allowed to remain between the stapler jaws until surgery.*

**Note:** *For optimal results the prepared PSD-V/stapler should be used within 60 minutes.*

17. Just prior to using the stapler, open the stapler and remove the foam spacer. (See Figure 4).

18. Visually inspect each buttress strip to ensure they are on the stapler jaws and cover the staple line.

**Note:** *Skip to step 27.*

### **LOADING STAPLER FOR ENDO GIA™ UNIVERSAL MODELS**

19. Position the open stapler onto the stapler loading unit (See Figure 2).

20. Ensure the anvil and cartridge sides of the loading unit are on the corresponding stapler jaws or the strips may not adhere to the stapler properly.

**Note:** *The cartridge and anvil sides of the PSD-V loading unit differ; substitution of one side for the other may interfere with alignment and adherence of the buttress strips.*

21. Close the stapler.

22. Apply light pressure to both sides of the stapler jaws to hold stapler into closed position.

23. Slide the assembly clip over the stapler jaws (See Figure 5).

24. Ensure that both stapler jaws are held firmly in place by the assembly clip and that the yellow guard, if present, is positioned between the assembly clip on the cartridge side.

**Note:** *PSD-V can be used immediately or allowed to remain between the stapler jaws until surgery.*

**Note:** *For optimal results the prepared PSD-V/stapler should be used within 60 minutes.*

25. Immediately before using, remove stapler loading unit (See Figure 6).

26. Visually inspect each buttress strip to ensure they are on the stapler jaws and cover the staple line.

### **IMPLANTING PSD-V INSTRUCTIONS FOR ALL STAPLER MODELS**

27. Follow *Instructions For Use* supplied by the stapler manufacturer.

28. If necessary, cut the end of the PSD-V to remove the dissected tissue (See Figure 7).

29. Discard any open PSD-V and Gel pouches. These cannot be re-sterilized or reused.

### **DISCLAIMER OF WARRANTIES:**

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. SSI will replace any device, which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Peri-Strips Dry is a registered trademark of Synovis Life Technologies, Inc., (US Patent #5,752,965; 6,312,474; 6,652,594; Additional Patents Pending)
A Division of Synovis Life Technologies, Inc.